

The Ocenco M-20.2 EEBD - THE WORLD'S SMALLEST AND MOST DURABLE EEBD.

Think of it as a Life Preserver You Wear on Your Belt


The Ocenco M-20.2 compressed oxygen EEBD (Emergency Escape Breathing Device) provides up to 32 minutes of protection.

The Ocenco M-20.2 EEBD can be donned in seconds; simply unlatch the case, pull out the unit, and insert the mouthpiece and nose clip.

The attached hood can be donned at anytime during the escape.

The compressed oxygen and mouthpiece combination allows the Ocenco M-20.2 EEBD to be donned in a smoke filled environment.

The compact Ocenco M-20.2 EEBD can be belt worn in all confined spaces.



The Ocenco M-20.2 EEBD shown donned and ready for escape.

The secondary container and optional mounting bracket shown in the foreground.

The M-20.2 Meets the New SOLAS Requirement for EEBD's

The Ocenco M-20.2 EEBD has been selected for exclusive use by the US Navy, US Coast Guard and the Hellenic Navy.







THE OCENCO M-20.2 EEBD:

Quick, easy use — pulling the unit from the case automatically starts oxygen flow.

Low life cycle cost — 15 year service life – annual shipboard visual inspection.

Belt wearable — so light and compact it can be worn comfortably on a belt.

Full visibility — the clear Teflon® hood protects the user from hazardous environments while allowing a full field of view.

Maximum Protection — Teflon* hood and breathing bag provide excellent heat and chemical resistance.

M-20.2 EEBD EMERGENCY ESCAPE BREATHING DEVICE Breathing bag Reliet valve M-20.2 PHYSICAL CHARACTERISTICS AND PERFORMANCE DATA



















| Typical performance duration | 15-32 minutes |
|--|---|
| Donning time | less than 10 seconds |
| Weight | 3.1 lbs (1.4 kg) |
| Dimensions | M-20.2 EEBD - $6.8 \times 6.8 \times 3$ in $(17 \times 17 \times 7.6$ cm) Secondary Container - $8.2 \times 9.2 \times 4.2$ in $(21 \times 23 \times 10$ cm) |
| Storage temperature range | -4°F to 149°F (-20° C to 65° C) |
| Service life | 15 Years |
| Oxygen delivery system | Automatic on, compressed oxygen, demand regulated |
| Inspection | Annual visual for stored units |
| M20.2 NSN No. M-20.2T Trainer NSN No. | NSN 1HM 4240-01-439-5937 NSN 1HM 4240-01-459-0078 |



LakeView Corporate Park 10225 82nd Avenue Pleasant Prairie, WI 53158-5801 U. S. A.

Phone: (262) 947-9000 Fax: (262) 947-9020

ACCESSORIES

M-20.2T Trainer — complete with two extra mouthpieces, in a light blue secondary container. Teaches donning from both the worn and stored position.

M-20.2 Secondary container — included with each M-20.2 but may be ordered separately.

M-20.2 Storage bracket — Bright orange welded steel bracket with quick opening Velcro[™] closure.